

Fralin room 004

Type of material:

Transport of material into room:

- Material must be placed inside primary container (tube/ plate, flask, etc.) in lab
- The tops to primary containers should be secured shut with parafilm or tape
- Disinfected on the outside
- Primary container placed in leakproof secondary container

Entering Room:

- You must bring your own gloves and any other supplies you need. All material must be handled with gloves on while in room.
- There will be 70% ethanol, Kimwipes, and hand sanitizer in the room
- Keep tops on primary containers unless it is absolutely necessary to remove them
- Disinfect (wipe with kimwipe saturated with ethanol) microscope stage when finished with procedure. Note: do not allow ethanol to get on the objective lenses.
- Gloves must be removed and you must use hand sanitizer before leaving the room
- Gloves must be taken back to your lab for disposal.
- Place all primary containers and gloves into leakproof secondary container

Exit room: